




AURORA HD

My dream to see the vibrant world again

Offering superior picture quality and a large viewing area, Aurora HD is ideal for reading larger books, writing or doing hobbies. You can also perform sorts of tasks easily with the foldable flexible swivel monitor in diverse directions. 24" widescreen brings more content.



-  Plug and play
-  Movable monitor
-  High definition image

Technical Specifications

Display	24" LCD widescreen, Resolution Ratio:1080 Full HD
Camera	High definition camera, Resolution ratio: 1280x720p
Magnification	2x to 70x
Color modes	17 modes
Default color modes	5 colors
Basic functions	Zoom in/out, Color mode adjustment, Brightness adjustment
Advanced functions	Freeze frame, Focus lock, Find feature, Reading lins and windows, OCR/ TTS function, Memory setting
Weight	Approximate 15.3 kg
Working height	18 cm
Certificates	CE FC RoHS

