

Butterfly 3.5"

For you with you

Pocket-size electronic magnifier,
Butterfly is designed for people on
the go.



Easy operation



5 color modes




Pocket-size on the go

Technical Specifications

| | |
|---------------|--|
| Display | 3.5" TFT screen |
| Magnification | 1.5X ~ 15X |
| Color modes | 5 modes |
| Image freeze | Yes |
| Battery | 3.5 hours Li-ion battery, rechargeable |
| Writing stand | Yes |
| Dimensions | 9.4cm×8cm×1.9cm |
| Weight | 120g |
| Accessories | AC charger, carry case, wrist strap |
| Certificates | CE; FCC; RoHS |



Contact us

 Tel: +86-571-87006308
Email: sales@zoomax.com
Web: www.zoomax.com