




# M5

## Miracle, at Your Fingertips!

M5, the first handheld video magnifier adopted **2HD** design - HD camera plus HD screen, in the global market, presents super crystal-clear image for all age group. Being proved by hundreds of testers from low vision organizations, schools and individuals, its freely touch and handy features allow everyone to experience the cutting-edge device you'll never imagine.



NEW

-  The **FIRST 2HD** design for 5" unit , HD CAMERA and HD SCREEN
-  The **LIGHTEST** weight 5" HD unit on the market
-  The **EASY** touch feature proved by testers globally

**ZOOMAX**<sup>®</sup>  
www.zoomax.com

# Technical Specifications

Display	5" TFT touch screen. Screen resolution: 1280x720p
HD Cameras	5 megapixels. Resolution: 1280x720p
Magnification	2.3X-16X
Color modes	10 contrast color + full color
Favorite color modes	Full color mode+ 2 favorite contrast color modes(reversed color)
Freeze/save photos	Yes, store up to 60 photos
Panning	Easy panning to view detailed parts of the photo
Reading stand	Yes, built-in reading stand
Distance viewing	Yes
Brightness adjustment	Yes, 6 brightness levels
Led light	Turn on/off
Connect to TV / monitor	Support HDMI
Continuous use	(Fully charged battery) approx. 4 hours battery life (depending on Mode and Brightness level settings).
Charging time	Approx. 4.5 hours
Weight & Dimensions	265g (9.35oz). 173mm* 80mm* 20mm(6.8"*3.1"*0.8")



## Contact us



Tel: +86-571-87006308  
 Email: sales@zoomax.com  
 Web: www.zoomax.com