

MARS **HD**

Explore new world

Perfect for students and professionals, this versatile video magnifier offers 3 viewing modes making it easy to view objects close-up, materials on your desk, or presentations across the room!



 3 viewing modes

 HD camera

 Less than 0.9kg


ZOOMAX[®]
www.zoomax.com

Technical Specifications

| | |
|----------------------|---|
| Camera | 12X optical zoom, autofocus HD with resolution 1280 x 720 pixels |
| Frame rate | 60 full digital frame /second @ 720p |
| Magnification | 2X ~ 80X (based on 14"screen) |
| Color modes | Full color + 16 high contrast color modes |
| 3 viewing modes | Close, distance and self-viewing |
| Ergonomic adjustment | Rotate 330° vertically or 180° horizontally for multi-angles reading |
| USB/HDMI connection | Connect to Mac, PC and monitor |
| Functions | Zoom in/out, Color mode adjustment, Brightness adjustment, Reading lines/windows, Freeze frame, Split screen, Find, Image saving, Led light on/off, Setting memory |
| Battery | Around 4 hours continuous use, rechargeable Li-ion battery |
| Dimensions (folded) | 22.6cm×8.6cm×12.5cm |
| Weight | Less than 0.9kg |
| Certificates | CE; FCC; RoHS |



Contact us



Tel: +86-571-87006308

Email: sales@zoomax.com

Web: www.zoomax.com