




# Panda **SD**

See Simpler, Live Easier

With a compact design and competitive price, Panda is perfect for a home or apartment with limited space and budget. Excellent features of Zoomax Panda allow the vast majority of visually impaired people to easily read books, do homework, and handle documents.



-  3 buttons only
-  Space-saving size
-  High quality with affordability

## Technical Specifications

Display	19" LCD screen
Camera	Auto focus
Magnification	2x to 70x
Color modes	17 modes
Favorite color mode	Full color + 4 favorite color modes
Default color modes	5 modes
Ergonomic adjustment	Screen height adjustable
Functions	Zoom in/out, Color mode adjustment, Brightness adjustment, Memory setting
Working height	25 cm
XY table dimensions	36.8cm x 36.6cm x 4.6cm X dimension travel: 20cm Y dimension travel: 23cm
Weight	17.6 kg
Certificates	CE FC RoHS

