


Snow 4.3"


The World is Clear Again

Zoomax Snow 4.3 inch handheld video magnifier is given advanced image technology and camera. It offers sharp image even at the lowest magnification and under fast movement. With large buttons and the unique one hand control design, it's very easy to use even without any learning at the very beginning. You can read and write on the go or at home.



 Easy to use with intuitive buttons

 Light weight

 Big and intuitive buttons

Technical Specifications

| | |
|-------------------|---|
| Display | 4.3" TFT screen |
| Camera | Auto focus |
| Magnification | 2x to 16x |
| Color modes | 10 modes |
| Image freeze | Yes |
| Writing stand | Yes |
| Connect to TV | PAL / NTSC |
| Battery | 3.5 hours Li-ion battery, replaceable & rechargeable |
| Power saving mode | Enter into power saving mode automatically when it has not been used for over 3 minutes |
| Weight | 204g (7.20oz) with battery / 165g (5.82oz) without battery |
| Dimensions | 160mm x 80mm x 24.8mm (6.30" x 3.15" x 0.98") |
| Certificates | CE FC RoHS |

