





Panda HD

Cheers for Yourself

With a high specification and the compact design, Zoomax Panda HD offers the most attractive advantages of a full line desktop electronic magnifiers, such as crisp image in high definition, intuitive use, space-saving, movable screen design etc.



-  Continuous crisp image at all level
-  Very easy to learn and use for anyone
-  Unexpected space-saving
-  High specification and better experience

Technical Specifications

Display	19" LCD screen
Camera	High defination, Auto focus, Resolution Ratio: 1280 x 720p
Magnification	2x to 70x
Color modes	17 modes
Favorite color mode	Full color + 4 favorite color modes
Default color modes	5 modes
Ergonomic adjustment	Screen height adjustable
Functions	Zoom in/out, Color mode adjustment, Brightness adjustment, Memory setting
Working height	26 cm
XY table dimensions	36.8cm x 36.6cm x 4.6cm X dimension travel: 20cm Y dimension travel: 23cm
Weight	Approximately 17.3 kg
Certificates	CE FC RoHS

