

Snow 10 Pro

HANDY SPEECH, TABLET-SIZE SCREEN,
FOR ACTIVE LIFESTYLE



"I use it to listen to recipes in the kitchen with my wife. Cooking is a fun activity again."

- Fred J. Henriksen

"I listen to books each evening before I go to sleep."

- David K. Tiedemann

"I read the newspaper again, by listening. I also take it with me on the train to my friend's home."

- Shawn L. Alonzo

THE ONLY 10-INCH VIDEO MAGNIFIER WITH OCR AND SPEECH

High quality, human-sounding speech, brings the printed word to life. Listen to print material you have scanned or with imported PDF files. Snow 10 Pro is lightweight and can be used virtually everywhere including restaurants, stores, while traveling, at school, at work or at the library.



IS THE TEXT-TO-SPEECH FUNCTION EASY TO LEARN?

"I am 69-years old and tried this at a conference. It took me just seconds to learn that reading is just a simple 2-step process." - Anna S. Marrero




HOW DOES SNOW 10 PRO COMPARE TO OTHER HANDHELD VIDEO MAGNIFIERS?

Besides being the first and only 10-inch magnifier with available OCR and speech, Snow 10 Pro is the only handheld video magnifier that includes voice prompts to help guide you while using the device.



WHY SNOW 10 PRO INSTEAD OF A DESKTOP VIDEO MAGNIFIER?

The lightweight design and included carry case make it convenient to carry around. Weighing just 2 pounds and running on battery for 3 hours on a full charge, mean you can use Snow 10 Pro virtually everywhere.



HANDY SPEECH



WHY ZOOMAX SNOW 10 PRO?

- High accuracy powered by world leading OCR supplier ABBYY
- Human-sounding speech brings the printed word to life
- User-friendly interfaces
- Unique voice prompt feature to guide users

IT READS FOR YOU

- Scan full-page printed materials and read aloud
- Read aloud imported PDF files
- Choose from 20 languages
- Start speech from the word you point at
- Adjust voice volume & speed

EVERYTHING IN A PACKAGE

- 13 megapixels camera
- 3 viewing modes including Distance & Selfie
- Full color & 10 high-contrast color modes
- Audio files playable
- Store up to 1,000 images
- Voice Prompt
- Weighing 2 pounds + 3 hours continuous use
- 2 hours fast charging

ACCESSORY

- Optional Zoomax easeScan with HD camera designed for OCR & Speech
- Protective carry case



ZOOMAX (USA) INC.

 TOLL-FREE: (866) 296-8388

 WWW.ZOOMAXUSA.COM

 SALES@ZOOMAXUSA.COM



LOCAL DEALER:

ZOOMAX[®]