

Distance Camera for Snow 12



A great companion for our portable video magnifier Snow 12, this optional distance camera is helpful in the classroom or in the workplace.



● **Plug and Play**

Connect the distance camera to Snow 12 and automatically switch to distance view mode.

● **Clear Display**

With a viewing distance up to 35 meters (almost 115 feet), objects or text on the other side of the room can be viewed.

● **Portable**

The distance camera is lightweight (335g/11oz) and easy to carry around.

Specifications:

- Camera: 5 megapixels
- Connection: USB Type-C
- Dimensions: 13.5cm×10cm×24.5cm/5.3"×3.9"×9.6"

☎ +86-571-8700-6308

🌐 www.zoomax.com

✉ sales@zoomax.com