


ZOOMAX®

# LUNA 8

Simple, useful and stylish



**NEW**



# Extreme simplicity for everyone

## BACK TO BASICS

Only 1 minute to learn, Luna 8 is an all-new BASIC BUT ESSENTIAL handheld electronic video magnifier from Zoomax. Using Luna 8 is simpler than using a calculator. Whether reading a book or newspaper at home, using it to study for school, or taking it with you to work, Luna 8 is the perfect companion.



### **Exceptional Image Quality**

No glare, no ghosting. With 13MP camera and full HD screen display, Luna 8 displays every word on the screen sharply.



### **Basic Functions Only**

Just good old-fashioned zoom in & out, high contrast to assist with reading, and a camera for taking photos to look at later.




### **Light and Portable**

500g (1.1lb) weight and extraordinary 4.5 hours battery life, the sleekly designed Luna 8 is your perfect portable magnifier.



# Technical Specifications

Magnification	2.5X - 19X
Screen	8-inch HD screen; resolution: 1280 x 800
Camera	auto focus; 13 megapixels
Color modes	full color and 10 high contrast color modes
Distance view	yes
Image storage	up to 200
Reading line & masks	yes
Battery life	more than 4.5 hours continuous use
Charging time	less than 3 hours
Weight	500g (1.1lbs)
Dimensions	218 x 143 x 11mm (8.6 x 5.6 x 0.4 in)

 +86-571-8700-6308

 [www.zoomax.com](http://www.zoomax.com)

 [sales@zoomax.com](mailto:sales@zoomax.com)

---