

LUNA S

Better and Happier



Portable and Lightweight



Crystal-clear Images



FM Radio in the Pocket



● Simple and Useful

Luna S is a 4.3-inch handheld electronic magnifier. It brings you crystal-clear images and FM radio feature in a compact body. Put Luna S in your pocket, and you are good to go: read menu, check bank statements, listen to radio, anything anywhere.

● Technical Specifications

- Magnification: 2x - 19x
- Screen: 4.3" TFT screen
- Color modes: 10 color modes
- Flashlight mode: Yes
- FM radio: Yes
- Battery life: 3h continuous use
- Foldable writing stand: Yes
- Weight: 7 oz (200g with battery / 167g without battery)
- Dimensions: 16cm x 8cm x 2.4cm /6.3" x 3.15" x 0.94"