


LUNA 8

Simple, useful and stylish





Extreme simplicity for everyone

8 Hours Battery Life

Only 1 minute to learn, Luna 8 is an all-new ESSENTIAL handheld electronic video magnifier from Zoomax. Using Luna 8 is simpler than using a calculator. Whether reading a book or newspaper at home, using it to study for school, or taking it with you to work, Luna 8 is the perfect companion. Its battery life is extremely long as 8 hours that enable you to get rid of charging the unit frequently.



Exceptional Image Quality

No glare, no ghosting. With 13MP camera and full HD screen display, Luna 8 displays every word on the screen sharply.




Easy-to-use Design

Luna 8 comes with zoom in & out, high contrast to assist with reading, and a camera for taking photos to look at later.



Light and Portable

500g (1.1lb) weight and extraordinary 8 hours battery life, the sleekly designed Luna 8 is your perfect portable magnifier.



Magnification	2.5X - 19X
Screen	full HD, 8-inch screen; screen resolution: 1920 x 1080
Camera	auto focus; 13 megapixels
Color modes	full color and 10 high contrast color modes
Distance view	yes
Image storage	up to 200
Reading line & masks	yes
Battery life	8 hours continuous use
Charging time	less than 3 hours
Weight	500g (1.1lbs)
Dimensions	218 x 143 x 11mm (8.6 x 5.6 x 0.4 in)

+86-571-8700-6308

www.zoomax.com

sales@zoomax.com

(USA) Toll-free: (833)265-3324

www.zoomaxusa.com

info@zoomax.com
