

LUNA HD 24 Pro

Enhance Your Vision Experience



- Plug and Play
- Flexible Movable Display
- User Friendly and Work Well
- HD Image, Accurate Color, and Brilliance





LUNA HD 24 Pro

The newly presented Luna HD 24 Pro is a stable, high-quality, and powerful full HD desktop video magnifier and provides the best solution for the legally blind and visually impaired.



Maximum Smooth in Reading

- Applies smooth algorithms to reduce the effects of aliasing.
- HD camera provides a wider field of view and superior image quality.
- Unique finding function makes it easy to keep track of location and prevent missing words or lines.



Maximum Comfort in Using

- 3 Physical Buttons for adjusting magnification, contrast, brightness.
- No installation. Unpackaging, plug it in, then use immediately.
- Simply put the materials on working surface, adjust the settings to your preference.



Maximum Flexible in Adjusting

- Flexible Rotating Experience.
 Humanized designs for various environments, always have the right reading angle and height.
- Available Working Area.
 Extra-large working surface,
 featuring easy-to-move platform with XY table.





Technical Specifications

Magnification	2.5x - 70x	
Color Mode	True color + 20 kinds of color mode	
Writing Height	Approx. 17.5cm	
Camera	Optical zoom lens, auto-focusing with resolution of 1280*720	
Screen	24-inch monitor with resolution of 1920*1080	
Weight	Approx. 15Kg	
Basic Function	Zoom in/ out, color mode adjustment, find function, contrast adjustment, brightness adjustment	
Advanced Function	20+ colors modes adjustment, find feature, reading line and mask, focus lock, freeze frame, memory setting	
Power Adapter	Input 100-240V, Output 19V/3A	
Certificates	CE F© RoHS	
+86-571-8700-6308 (833)265-3324 (USA)	www.zoomax.com www.zoomaxusa.com (USA)	sales@zoomax.com info@zoomaxusa.com (USA)