

Acesight

Wear, See, and Enjoy

Electronic glasses help visually impaired people regain visual independence



















Acesight VR is lightweight, offers superior image quality, and is easy to operate thanks to the wireless 3-button controller. It's an excellent choice for people with low vision.

- Easy to use: easy operation via a wireless controller
- Incredible image quality: 48 megapixel camera and full HD resolution
- No image lag: best-in-class 60 FPS frame rate
- Comfy to wear: lightweight and cool operating temperature

Acesight & Acesight S

Acesight and Acesight S are FDA approved Class I medical devices that deliver significant improvement in near vision ability for people with low vision.

Acesight & Acesight S Comparison

Function	Acesight	Acesight S
Zoom in/out		
High Contrast		
Outline		
Controller Lock		
Focus Lock		
Floating Reading Mode		
Narrow Mode		—
Find Function		—



Mobility

The open design means you can walk while wearing Acesight.



Options


RP or limited FOV – Acesight.
For AMD – Acesight S.



Long Use

Great battery life provides for 4.5 hours continuous use.

WHO CAN BENEFIT FROM ACESIGHT FAMILY?



The Acesight family of head mounted displays were designed to deliver greater visual independence for people with a visual acuity between 20/100 and 20/600.

Suitable for the following eye conditions:

- ✓ Macular Degeneration
- ✓ Glaucoma
- ✓ Diabetic Retinopathy
- ✓ Stargardt's Disease
- ✓ Nystagmus
- ✓ Esotropia
- ✓ Retinitis Pigmentosa
- ✓ And more



Perfect for watching TV



Computer screens look great



Perform daily tasks

Please contact our sales team for more information.

 www.acesight.com

 info@acesight.com